

Anatomy: **ABDOMEN**

- Exams

Sub-Anatomy: **Eovist Liver BH**

- Routine

Coil: **SENSE XL TORSO**

SEQUENCE - BASICS				VARIABLE PARAMETERS														
PLANE	SEQ	Slice / Gap (mm)	Misc / Comment	MTX	FOV (cm)	% RFOV	Scan %	Pixel size (mm)	TR	TE	FA/TI	Nex Avg Acq	NS	ETL Turbo Factor	Phase Encode	Scan TIME (min)	Pixel Shift BW-kHz	
<b>ROUTINE</b>																		
COR	2D T2 SSFSE	5/1	BH															
AX	2D T1 IP/OP	5/1	BH															
AX	3D mDIXON DYN	4 -> 2	Pre, Art,40,90,120															
AX	2D T2 SSFSE	5/1	BH															
AX	2D T2 FS SSFSE	5/1	BH															
AX	2D DWI FS	7/1	b=0,50,400,800															
AXTOP	3D mDIXON 20MIN	3 -> 1.5	Single phase Cover top half of liver Allow table movement					1.3x1.63 -> 0.78x0.77									15.1 sec	
AXBOT	3D mDIXON 20MIN	3 -> 1.5	Single phase Cover bottom half of liver Allow table movement															
<b>↓ OPTIONAL ↓</b>																		

**ORDERABLE: MR ABDOMEN WWO CONTRAST**

TIPS: (1) Cover whole liver.

(2) Make sure the "3D mDIXON 20MIN" is scanned at least 20 min after Eovist injection (please wait if you have to).

(3) Make sure patient has a good breath-hold for the "3D mDIXON 20MIN" scans. Use three breathing instructions if needed.

(4) When planning last two scans, make sure "AXTOP" covers top-half of the liver, "AXBOT" covers bottom-half, with some slice overlapping between the two scans.