

Anatomy: **Pelvis**
 Sub-Anatomy: **Pelvis 1.5T or 3T**

- Exams **ORDERABLE- Pelvis**
 - Routine Coil: **Torso coil**


SEQUENCE - BASICS																	
PLANE	SEQ	Slice thickness (mm)	Misc / Comment	M T X	% R F O V	Gap (mm)	Voxel size (mm)	TR	TE	F A/ TI	Nex Avg Acq	NS	ETL Turbo Factor	Phase Encode	Scan TIME (min)	Pixel Shift BW-kHz	
ROUTINE																	
	3 plane scout		Only GRE														
1	AX T2 FS	4x0.6x0.7				0.4		4000	50-65								
2	Ax T1	4x0.6x0.7				0.4		600	6-9								
3	Cor T1	5x0.6x0.7				0.5		600	6-9								
4	Cor STIR	5x0.6x0.7				0.5		2000	25-35								
5	Sagittal T2 Dixon	5x0.6x0.7				0.5		4000	50-65								
↓ OPTIONAL ↓																	
	Ax STIR	4x0.6x0.7	Failed fat sat			0.4		2000	25-35								

Instructions: FOV and Coverage- On axials and coronals, cover from L4-5 to just below lesser trochanters and from skin to skin. On sagittals, cover from greater trochanter to greater trochanter.


Others- Always put marker at the site of most pain, injury or swelling. If motion observed, stop the scan, instruct patient to lay still and repeat the sequence.



Axial Coverage



Coronal Coverage



Sagittal Coverage