

Anatomy: UPPER EXTREMITY  
 Sub-Anatomy: Wrist- 1.5T Ortho

- Exams **ORDERABLE- WRIST**  
 - Routine Coil: Wrist/Flex coil


SEQUENCE - BASICS																
PLANE	SEQ	Slice thickness (mm)	Misc / Comment	MT X	Gap	Slices	Sc an %	TR	TE			NS	ETL Turbo Factor	Phase Encode	Scan TIME (min)	Pixel Shift BW-kHz
<b>ROUTINE</b>																
	3 plane scout		Only use GRE													
1	Cor PD	3x0.4x0.5			10%			4000	40-45							
2	Cor PD FS	3x0.4x0.5			10%			4000	35-40							
3	Ax PD	3x0.4x0.5			10%			4000	40-45							
4	Ax PD FS	3x0.4x0.5			10%			4000	35-40							
5	Sag PD FS	3x0.4x0.5			10%			4000	35-40							
6	Cor T2 Dixon	4x0.5x0.6			10%											
<b>↓ OPTIONAL ↓</b>																
	AX STIR	4x0.4x0.5	Failed fat sat		10%			4000	30-35							

**Instructions: FOV and Coverage-** On coronal, cover from skin to skin. On axials, cover radial meta-diaphyseal junction to bases of metacarpals. On sagittals, cover from skin to skin. Superman position or wrist at side if cannot tolerate supermanposition


**Others-** Always put marker at the site of most pain, injury or swelling. If motion observed, stop the scan, instruct patient to lay still and repeat the sequence.



Axial



Coronal



Sagittal